MPFD Foil Activation Experiment Resource

John Boyington Kansas State University

For each wand:


Inventory

- 1. Loaded wand
- 2. Stopwatch
- 3. Gloves
- 4. Sample Bags (Four plastic bags labeled #1, #2, #3, and #4; #1 corresponds to the highest axial position.)
- 5. Sample pig for transportation
- 6. Neutron dose monitor.

Procedures

- 1. Raise reactor power to P then insert wand into reactor core.
- 2. Irradiate foil for t_i .
- 3. Scram reactor and promptly remove wand from reactor core.
- 4. While dose is above 50 mR/h at 1 ft., keep stored in fuel rod storage.
- 5. Pull wand to surface and remove internals. Remove foils from internals and place in labeled bags. Place bags in sample pig.
- 6. Count foils one at a time with HPGe for listed counting time.

Principle Reactions




Reaction	$T_{1/2}$	ROI (eV)	Important Gammas (keV)
$^{115}\mathrm{In}(\mathrm{n},\gamma)^{116}\mathrm{In}$	54 m	7.0021e-03, 1.6130e+00	417, 819, 1090, 1293, 1508, 2111
115 In(n, γ) 116 In Cd	54 m	1.1955e+00, 1.9916e+00	417, 819, 1090, 1293, 1508, 2111
$^{197}\mathrm{Au}(\mathrm{n},\gamma)^{198}\mathrm{Au}$	2.7 d	6.7266e-03, 5.2684e+00	412, 676, 1088
197 Au(n, γ) 198 Au Cd	2.7 d	4.0752e+00, 7.1730e+00	412, 676, 1088
$^{103}\text{Rh}(\text{n,n'})^{103m}\text{Rh}$	56.12 m	4.4469e+05, 5.1947e+06	40
$^{27}\mathrm{Al}(\mathrm{n},\alpha)^{24}\mathrm{Na}$	15.03 h	6.4564e+06, 1.1695e+07	1369, 2754

Indium

Power Level: 1.0 kW(th) Time at Power: 60 s Wait Time: 400 s

Total Activity at Removal: 1.79e+02 μCi

Position	Mass mg	Start Counting s	Counting Time s	Counting Activity μCi
1	1.7	460	300	4.51e+01
2	1.5	760	300	3.73e+01
3	1.4	1060	300	3.27e+01
4	1.6	1360	300	3.50e+01




Reaction	$T_{1/2}$	ROI (eV)	Important Gammas (keV)
(n,γ)	54.0 m	7.00e-03, 1.61e+00	138(0.03), 417(0.36), 819(0.17), 1090(0.53), 1293(0.8), 1508(0.11),
(n, n')	4.4 h	1.33e+06, 5.96e+06	335(0.5)

Indium (Cd)

Power Level: 1.0 kW(th) Time at Power: 60 s Wait Time: 400 s

Total Activity at Removal: 8.85e+01 μCi

Position	Mass mg	Start Counting s	Counting Time s	Counting Activity μCi
1	1.7	460	300	2.23e+01
2	1.5	760	300	1.84e+01
3	1.4	1060	300	1.61e+01
4	1.6	1360	300	1.73e + 01




Reaction	$T_{1/2}$	ROI (eV)	Important Gammas (keV)
(n, γ)	54.0 m	1.20e+00, 1.99e+00	138(0.03), 417(0.36), 819(0.17), 1090(0.53), 1293(0.8), 1508(0.11),
(n, n')	4.4 h	1.34e+06, 5.97e+06	335(0.5)

Gold

Power Level: 100 kW(th) Time at Power: 45 s Wait Time: 300 s

Total Activity at Removal: 2.73e+03 μCi

Position	Mass mg	Start Counting s	Counting Time s	Counting Activity μCi
1	5.0	345	300	4.24e+01
2	4.35	645	300	3.69e+01
3	4.3	945	300	3.64e+01
4	4.37	1245	300	3.70e+01




Reaction	$T_{1/2}$	ROI (eV)	Important Gammas (keV)
(n,γ)	2.7 d	6.73e-03, 5.27e+00	412(0.95), 676(0.01), 1088(0.002)
(n,2n)	6.2 d	8.83e+06, 1.35e+07	333(0.25), 356(0.94), 426(0.06), 1091(0.002)
(n, n')	7.8 s	3.24e+05, 5.25e+06	130(0.08), 279(0.75)

Gold (Cd)

Power Level: 100 kW(th) Time at Power: 45 s Wait Time: 300 s

Total Activity at Removal: 2.49e+03 μCi

Position	Mass mg	Start Counting s	Counting Time s	Counting Activity μCi
1	5.0	345	300	2.28e+01
2	4.35	645	300	1.98e+01
3	4.3	945	300	1.96e + 01
4	4.37	1245	300	1.99e + 01




	Reaction	$T_{1/2}$	ROI (eV)	Important Gammas (keV)
ſ	(n, γ)	2.7 d	4.08e+00, 7.17e+00	412(0.95), 676(0.01), 1088(0.002)
ſ	(n,2n)	6.2 d	8.83e+06, 1.35e+07	333(0.25), 356(0.94), 426(0.06), 1091(0.002)
	(n, n')	7.8 s	3.28e+05, 5.28e+06	130(0.08), 279(0.75)

Rhodium

Power Level: 100 kW(th) Time at Power: 600 s Wait Time: 5000 s

Total Activity at Removal: 2.20e+05 μCi

Position	Mass mg	Start Counting s	Counting Time s	Counting Activity μCi
1	0.7	5600	600	7.13e+00
2	0.55	6200	600	4.88e+00
3	0.5	6800	600	3.91e+00
4	0.55	7400	600	3.79e + 00




Reaction	$T_{1/2}$	ROI (eV)	Important Gammas (keV)
(n, γ)	4.4 m	4.67e-03, 1.39e+00	51(0.47), 78(0.025), 560(0.026), 770(0.0018)
(n, n')	56.1 m	4.45e+05, 5.19e+06	40(0.004)

Aluminum

Power Level: 250 kW(th) Time at Power: 3600 s Wait Time: 1800 s

Total Activity at Removal: 9.62e+02 μCi

Position	Mass mg	Start Counting s	Counting Time s	Counting Activity μCi
1	0.3	5400	1800	1.92e-01
2	0.2	7200	1800	1.75e-02
3	0.1	9000	1800	3.83e-03
4	0.2	10800	1800	6.44e-03



	Reaction	$T_{1/2}$	ROI (eV)	Important Gammas (keV)
	(n,p)	9.5 m	3.42e+06, 9.14e+06	180(0.007), 840(0.7), 1013(0.3)
Г	(n, α)	15.0 h	6.46e+06, 1.17e+07	1369(1), 2754(1)
	(n, γ)	2.2 m	2.82e-03, 4.15e-01	1780(1)

Useful Links

Activation Calculator

https://www.ncnr.nist.gov/resources/activation/

Online Spectrum Catalogs for Ge and Si(Li)

http://www4vip.inl.gov/gammaray/catalogs/ge/catalog_ge.shtml

Decay Radiation Search

https://www.nndc.bnl.gov/nudat2/indx_dec.jsp

Evaluated Nuclear Data File (ENDF) Retrieval & Plotting https://www.nndc.bnl.gov/sigma/